

Extend 180™ Long-Term Absorbable Synthetic Implants


6 Months of Dry Eye Relief*

Extend 180 Long-Term Dissolvable Implants:

- Provide occlusion for 6 full months of effective dry eye relief
- Are made of polydioxanone material
- Are available in 3 sizes and 2 packaging configurations (single pair and 10-pair box)

Indications:

- Post-ocular surgery
- Seasonal dry eye
- Contact lens intolerance
- Dry eye associated with digital eye strain



For ordering information call:
866-906-8080 customersupport@beaver-visitec.com

Extend 180 Absorbable Synthetic Implants			
	Description	Item #	Quantity
	0.3 mm X 2.0 mm Single Pair	A13-2033SP	1 pair box, 2 occluders total
	0.4 mm X 2.0 mm Single Pair	A13-2044SP	
	0.5 mm X 2.0 mm Single Pair	A13-2055SP	
	0.3 mm X 2.0 mm Bulk Pack	A13-2033	10 pair box, 20 occluders total
	0.4 mm X 2.0 mm Bulk Pack	A13-2044	
	0.5 mm X 2.0 mm Bulk Pack	A13-2055	
	Extend 180 Absorbable Synthetic Implant Variety Pack	A13-2000	10 pair box of: 2 pairs — 0.3 mm X 2.0 mm 5 pairs — 0.4 mm X 2.0 mm 3 pairs — 0.5 mm X 2.0 mm

Beaver-Visitec International, Inc.
411 Waverley Oaks Road, Waltham, MA 02452
Tel: 1.866.906.8080 Fax: 1.866.906.4304
www.beaver-visitec.com



* 510(k) Summary K162361

BVI, BVI Logo and all other trademarks (unless noted otherwise) are property of Beaver-Visitec International ("BVI") © 2017 BVI.

590519-06-17